

## Letter of Confirmation

Topic : Covid Variants Detection including Omicron variant(B.1.1.529)  
by the Sugentech Rapid Ag Test Kit (SGTi-flex COVID-19 Ag)

To. whom it may concern,

We, Sugentech, the manufacturer of the 'SGTi-flex COVID-19 Ag' would like to confirm the followings.

As our product, SGTi-flex COVID-19 Ag, only detects the SARS-CoV-2 "nucleocapsid protein", the product performance is not affected by variants with spike protein mutations.

Our product can detect the recombinant nucleocapsid protein of the SARS-CoV-2 variants such as Alpha (B.1.1.7), Beta (B.1.351), Gamma (P.1), Kappa (B.1.617.1), B.1.617, Delta (B.1.617.2), B.1.617.3, Epsilon (B.1.427 and B.1.429), and the new variant Omicron (B.1.1.529).

We hereby certify that the above information is true and correct.  
Please contact us if you have any further questions.

Sincerely yours,

November 30, 2021



Mijin Sohn/ CEO & President



SUGENTECH, INC.